


DLC1 Antibody (R36094)

Catalog No.	Formulation	Size
R36094-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10395
Applications	Western Blot : 1-5ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This DLC1 antibody is available for research use only.



IHC testing of FFPE human adrenal gland with DLC1 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Deleted in liver cancer 1, ARHGAP7; DLC1

Application Notes

Optimal dilution of the DLC1 antibody should be determined by the researcher.

1. This DLC1 antibody is specific to human isoforms 1 and 3, and will not detect isoform 2.

Immunogen

Amino acids SVAIRKRSWEEHC were used as the immunogen for this DLC1 antibody.

Storage

Aliquot and store the DLC1 antibody at -20oC.