


DLC1 Antibody (R36094)

Catalog No.	Formulation	Size
R36094-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10395
Applications	Western Blot : 1-5ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This DLC1 antibody is available for research use only.



IHC testing of FFPE human adrenal gland with DLC1 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Deleted in liver cancer 1, ARHGAP7; DLC1

Application Notes

Optimal dilution of the DLC1 antibody should be determined by the researcher.

1. This DLC1 antibody is specific to human isoforms 1 and 3, and will not detect isoform 2.

Immunogen

Amino acids SVAIRKRSWEEHC were used as the immunogen for this DLC1 antibody.

Storage

Aliquot and store the DLC1 antibody at -20oC.