


DLL1 Antibody (R33191)

Catalog No.	Formulation	Size
R33191-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	28514
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This DLL1 antibody is available for research use only.



Western blot testing of mouse heart lysate with DLL1 antibody at 0.3ug/ml. Predicted molecular weight ~78 kDa.

Description

Additional name(s) for this target protein: Delta-like 1

Application Notes

Optimal dilution of the DLL1 antibody should be determined by the researcher.

Immunogen

Amino acids ATQRHLTVGEEWSQD were used as the immunogen for this DLL1 antibody.

Storage

Aliquot and store the DLL1 antibody at -20oC.