


## DLL4 Antibody (Center Region) (R34555)

Catalog No.	Formulation	Size
R34555-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54567
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This DLL4 antibody is available for research use only.



Western blot testing of human [A], mouse [B] and rat [C] lung lysates with DLL4 antibody at 0.3ug/ml. Predicted/observed molecular weight 75~85kDa.

## Description

Additional name(s) for this target protein: Delta ligand 4

## Application Notes

Optimal dilution of the DLL4 antibody should be determined by the researcher.

## Immunogen

Amino acids KNTNQKKELEVDC were used as the immunogen for this DLL4 antibody.

## Storage

Aliquot and store the DLL4 antibody at -20oC.