


## DMRT2 Antibody (R34631)

Catalog No.	Formulation	Size
R34631-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10655
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This DMRT2 antibody is available for research use only.



IHC testing of FFPE human kidney with DMRT2 antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

### Description

Additional name(s) for this target protein: Doublesex and mab-3 related transcription factor 2

### Application Notes

Optimal dilution of the DMRT2 antibody should be determined by the researcher.

1. This DMRT2 antibody is specific to isoform 2.

## Immunogen

Amino acids KRPSSAITRVSQ were used as the immunogen for this DMRT2 antibody.

## Storage

Aliquot and store the DMRT2 antibody at -20oC.