


## DNAJB9 Antibody (R34965)

Catalog No.	Formulation	Size
R34965-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4189
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This DNAJB9 antibody is available for research use only.



Western blot testing of human liver lysate with DNAJB9 antibody at 1ug/ml. Predicted molecular weight: ~26kDa.

## Description

Additional name(s) for this target protein: DnaJ (Hsp40) homolog, subfamily B, member 9

## Application Notes

Optimal dilution of the DNAJB9 antibody should be determined by the researcher.

## Immunogen

Amino acids PDAEAKFREIAEAYE were used as the immunogen for this DNAJB9 antibody.

## Storage

Aliquot and store the DNAJB9 antibody at -20oC.