


DNAJB9 Antibody (R34965)

Catalog No.	Formulation	Size
R34965-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4189
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This DNAJB9 antibody is available for research use only.



Description

Additional name(s) for this target protein: DnaJ (Hsp40) homolog, subfamily B, member 9

Application Notes

Optimal dilution of the DNAJB9 antibody should be determined by the researcher.

Immunogen

Amino acids PDAEAKFREIAEAYE were used as the immunogen for this DNAJB9 antibody. **Storage** Aliquot and store the DNAJB9 antibody at -20oC.