



Dnajc5 Antibody (R35843)

Catalog No.	Formulation	Size
R35843-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	80331
Localization	Cytoplasmic
Applications	Western Blot : 0.5-1.5ug/ml IHC (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Dnajc5 antibody is available for research use only.





IHC testing of FFPE human skeletal muscle tissue with Dnajc5 antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: DnaJ homolog subfamily C member 5

Application Notes

Optimal dilution of the Dnajc5 antibody should be determined by the researcher.

Immunogen

Amino acids SATETTQLTADSHPS were used as the immunogen for this Dnajc5 antibody.

Storage

Aliquot and store the Dnajc5 antibody at -20oC.