


## Dnajc5 Antibody (R36261)

Catalog No.	Formulation	Size
R36261-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	80331
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Dnajc5 antibody is available for research use only.



Western blot testing of human cerebellum lysate with Dnajc5 antibody at 0.1ug/ml.  
Predicted molecular weight: 22 kDa.

## Description

Additional name(s) for this target protein: DnaJ (Hsp40) homolog, subfamily C, member 5

## Application Notes

Optimal dilution of the Dnajc5 antibody should be determined by the researcher.

## Immunogen

Amino acids TQLTADSHPSYHTD were used as the immunogen for this Dnajc5 antibody.

## **Storage**

Aliquot and store the Dnajc5 antibody at -20oC.