


# **DNAJC5G Antibody (R33631)**

Catalog No.	Formulation	Size
R33631-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	285126
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This DNAJC5G antibody is available for research use only.



IHC testing of FFPE human testis with DNAJC5G antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

# **Description**

Additional name(s) for this target protein: DNAJ (Hsp40) homolog, subfamily C, member 5 gamma

## **Application Notes**

Optimal dilution of the DNAJC5G antibody should be determined by the researcher.

### **Immunogen**

Amino acids DSGRKYQQNVQSQP were used as the immunogen for this DNAJC5G antibody.

# **Storage** Aliquot and store the DNAJC5G antibody at -20oC.