


DNAJC5G Antibody (R33631)

Catalog No.	Formulation	Size
R33631-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	285126
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This DNAJC5G antibody is available for research use only.



IHC testing of FFPE human testis with DNAJC5G antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: DNAJ (Hsp40) homolog, subfamily C, member 5 gamma

Application Notes

Optimal dilution of the DNAJC5G antibody should be determined by the researcher.

Immunogen

Amino acids DSGRKYQQNVQSQP were used as the immunogen for this DNAJC5G antibody.

Storage

Aliquot and store the DNAJC5G antibody at -20oC.