


DOK3 Antibody / Docking protein 3 (R36264)

Catalog No.	Formulation	Size
R36264-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	79930
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This DOK3 antibody is available for research use only.



Western blot testing of human PBMC lysate with DOK3 antibody at 0.1 ug/ml. Predicted molecular weight: ~53kDa, observed here at ~60kDa.

Description

Additional name(s) for this target protein: Docking Protein 3

Application Notes

Optimal dilution of the DOK3 antibody should be determined by the researcher.

Immunogen

Amino acids TRERRKGPAPCDRP were used as the immunogen for this DOK3 antibody.

Storage

Aliquot and store the DOK3 antibody at -20oC.