


## DPM3 Antibody (R34166)

Catalog No.	Formulation	Size
R34166-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q9P2X0
Gene ID	54344
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This DPM3 antibody is available for research use only.



Western blot testing of K562 lysate with DPM3 antibody at 2ug/ml. Predicted molecular weight: ~10/13kDa (short/long isoform).

## Description

Additional name(s) for this target protein: Dolichyl-phosphate mannosyltransferase polypeptide 3

## Application Notes

Optimal dilution of the DPM3 antibody should be determined by the researcher.

## **Immunogen**

Amino acids HDCEDAARELQSQ were used as the immunogen for this DPM3 antibody. This sequence is common to the short and long isoform.

## **Storage**

Aliquot and store the DPM3 antibody at -20oC.