


Drosophila E2F1 Antibody (R36191)

Catalog No.	Formulation	Size
R36191-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Drosophila
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q27368
Applications	Immunofluorescence : suitable ELISA (peptide) LOD : 1:32000
Limitations	This Drosophila E2F1 antibody is available for research use only.



Drosophila E2F1 antibody staining of the Drosophila Eye-Antennal Imaginal Disc (0.5ug/ml).

Description

Additional name(s) for this target protein: E2F transcription factor, drosE2F1

Application Notes

Optimal dilution of the Drosophila E2F1 antibody should be determined by the researcher.

Immunogen

Amino acids TKGLFQTPYRSARN were used as the immunogen for this Drosophila E2F1 antibody.

Storage

Aliquot and store the Drosophila E2F1 antibody at -20oC.