


DSCAM Antibody (R33651)

Catalog No.	Formulation	Size
R33651-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1826
Applications	Western Blot : 2-6ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This DSCAM antibody is available for research use only.



Western blot testing of human cerebellum lysate with DSCAM antibody at 2ug/ml.
Predicted molecular weight: ~174 kDa.

Description

Additional name(s) for this target protein: Down syndrome cell adhesion molecule

Application Notes

Optimal dilution of the DSCAM antibody should be determined by the researcher.

Immunogen

Amino acids DSWDSAQRTKDVSPQ were used as the immunogen for this DSCAM antibody.

Storage

Aliquot and store the DSCAM antibody at -20oC.