

DSG1 Antibody / Desmoglein 1 [clone MSVA-544M] (V6068)

Catalog No.	Formulation	Size
V6068-100UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V6068-20UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG2a, kappa
Clone Name	MSVA-544M
UniProt	Q02413
Localization	Cell junction, Cell membrane, Desmosome
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This DSG1/Desmoglein 1 antibody is available for research use only.

Application Notes

1. Optimal dilution of the DSG1/Desmoglein 1 antibody should be determined by the researcher.
2. This DSG1/Desmoglein 1 antibody is recombinantly produced by expression in CHO cells.

Immunogen

Bovine desmoglein (desmosomal glyco protein I) from nasal epidermis was used as the immunogen for the DSG1/Desmoglein 1 antibody.

Storage

DSG1/Desmoglein 1 antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

