

DSG3 Antibody / Desmoglein 3 [clone MSVA-543M] (V6069)

Catalog No.	Formulation	Size
V6069-100UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V6069-20UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, kappa
Clone Name	MSVA-543M
UniProt	P32926
Localization	Cell junction, Cell membrane, Desmosome
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This DSG3/Desmoglein 3 antibody is available for research use only.

Application Notes

1. Optimal dilution of the DSG3/Desmoglein 3 antibody should be determined by the researcher.
2. This DSG3/Desmoglein 3 antibody is recombinantly produced by expression in human HEK293 cells.

Immunogen

A recombinant fragment (around amino acids 379-491) human DSG3 protein (exact sequence is proprietary) was used as the immunogen for the DSG3/Desmoglein 3 antibody.

Storage

DSG3/Desmoglein 3 antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

