


## DTX3 Antibody (R34548)


Catalog No.	Formulation	Size
R34548-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	196403
Localization	Cytoplasmic
Applications	Western Blot : 0.3-1ug/ml IHC (FFPE) : 10ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This DTX3 antibody is available for research use only.



Western blot testing of NIH3T3 lysate with DTX3 antibody at 0.3ug/ml. Predicted molecular weight: ~38 kDa, observed here at ~90 kDa.



IHC testing of FFPE human brain cortex with DTX3 antibody at 10ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

## Description

Additional name(s) for this target protein: Deltex 3

## Application Notes

Optimal dilution of the DTX3 antibody should be determined by the researcher.

## Immunogen

Amino acids KNKVTVSKPVWD were used as the immunogen for this DTX3 antibody.

## Storage

Aliquot and store the DTX3 antibody at -20oC.