


## DUT Antibody (R35363)

Catalog No.	Formulation	Size
R35363-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1854
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This DUT antibody is available for research use only.



Western blot testing of HEK293 [A] and HeLa [B] lysates with DUT antibody at 0.1ug/ml.  
Predicted molecular weight: 27/18kDa (isoforms 1/2).

## Description

Additional name(s) for this target protein: Deoxyuridine triphosphatase

## Application Notes

Optimal dilution of the DUT antibody should be determined by the researcher.

## Immunogen

Amino acids QLRFARLSEHATAPT were used as the immunogen for this DUT antibody. This sequence is common to isoforms 1 and 2.

## Storage

Aliquot and store the DUT antibody at -20oC.