


DYDC1 Antibody (R35688)

Catalog No.	Formulation	Size
R35688-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q8WWB3
Gene ID	143241
Applications	Western Blot : 1-5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This DYDC1 antibody is available for research use only.



Western blot of HEK293 lysate overexpressing human DYDC1-MYC lysate (Lane A) and mock-transfected lysate (Lane B) with DYDC1 antibody at 1ug/ml. Lane C - same transfected lysate tested with anti-MYC tag (1:1000). Predicted molecular weight: ~21/22kDa (isoforms 1/2).

Description

Additional name(s) for this target protein: DPY30 domain containing 1

Application Notes

Optimal dilution of the DYDC1 antibody should be determined by the researcher.

Immunogen

Amino acids RNEDILHSEEATLD were used as the immunogen for this DYDC1 antibody. **Storage** Aliquot and store the DYDC1 antibody at -20oC.