


DYDC1 Antibody (R35688)

Catalog No.	Formulation	Size
R35688-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q8WWB3
Gene ID	143241
Applications	Western Blot : 1-5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This DYDC1 antibody is available for research use only.



Western blot of HEK293 lysate overexpressing human DYDC1-MYC lysate (Lane A) and mock-transfected lysate (Lane B) with DYDC1 antibody at 1ug/ml. Lane C - same transfected lysate tested with anti-MYC tag (1:1000). Predicted molecular weight: ~21/22kDa (isoforms 1/2).

Description

Additional name(s) for this target protein: DPY30 domain containing 1

Application Notes

Optimal dilution of the DYDC1 antibody should be determined by the researcher.

Immunogen

Amino acids RNEDILHSEEATLD were used as the immunogen for this DYDC1 antibody.

Storage

Aliquot and store the DYDC1 antibody at -20oC.