

DZIP1 Antibody / Cilium assembly protein DZIP1 Antibody (FY13601)

Catalog No.	Formulation	Size
FY13601	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Buffer	Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ .
UniProt	Q86YF9
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : 2-5ug/ml ELISA : 0.1-0.5 ug/ml
Limitations	For research use only, this DZIP1 Antibody / Cilium assembly protein DZIP1 Antibody is not intended for diagnostic or therapeutic use.

DZIP1 Antibody by IHC in rat testis. Immunohistochemistry of paraffin-embedded rat testis using Cilium assembly protein DZIP1 Antibody reveals brown staining in the seminiferous tubules, including germinal epithelial cells. Sections underwent heat-mediated antigen retrieval in EDTA buffer (pH 8.0), blocking with goat serum, and primary antibody incubation at 2 ug/mL overnight at 40C. Bound antibody was detected using an HRP-conjugated secondary antibody and DAB chromogen.

DZIP1 Antibody by WB in human and rodent samples. Western blot analysis using Cilium assembly protein DZIP1 Antibody detects a prominent immunoreactive band near 113 kDa in human HEL (lane 1), SH-SY5Y (lane 2), 293T (lane 3), rat testis (lane 4), and mouse testis (lane 5). The observed migration is compared with the calculated molecular weight of approximately 99 kDa. Reduced lysates containing 30 ug protein per lane were separated by 8% SDS-PAGE and transferred to nitrocellulose. After blocking with 5% nonfat milk, the membrane was incubated with primary antibody at 0.5 ug/mL overnight at 40C, followed by HRP-conjugated secondary antibody and chemiluminescent detection.

DZIP1 Antibody by IHC in human colon carcinoma. Immunohistochemistry of paraffin-embedded human colon carcinoma using Cilium assembly protein DZIP1 Antibody reveals heterogeneous brown cytoplasmic and perinuclear staining in malignant glandular epithelium. Sections underwent heat-mediated antigen retrieval in EDTA buffer (pH 8.0), blocking with goat serum, and primary antibody incubation at 2 ug/mL overnight at 4oC. Bound antibody was detected using an HRP-conjugated secondary antibody and DAB chromogen.

DZIP1 Antibody by IHC in human testicular carcinoma. Immunohistochemistry of paraffin-embedded human testicular carcinoma using Cilium assembly protein DZIP1 Antibody reveals brown immunoreactivity in scattered tumor cells. Sections underwent heat-mediated antigen retrieval in EDTA buffer (pH 8.0), blocking with goat serum, and primary antibody incubation at 2 ug/mL overnight at 4oC. Bound antibody was detected using an HRP-conjugated secondary antibody and DAB chromogen.

Description

DZIP1 Antibody / Cilium assembly protein DZIP1 Antibody detects DZIP1, also called cilium assembly protein DZIP1, a protein involved in the formation and function of primary cilia. Primary cilia are microtubule-based projections that organize signaling pathways at the cell surface. Their assembly requires coordinated activity at the mother centriole, where a basal body anchors the developing axoneme. DZIP1 has been implicated in events linking centriolar structures to the recruitment and organization of proteins needed for ciliogenesis.

Ciliary assembly depends on membrane docking, distal appendage function, and regulated trafficking of proteins into the ciliary compartment. DZIP1 has been studied alongside other centrosomal and ciliary regulators that help establish the conditions for axoneme growth. Its activity is particularly relevant to the formation of a signaling-competent cilium, not simply to the presence of microtubules. Defects in early assembly steps can disrupt the localization of receptors and signaling proteins that normally concentrate within the ciliary membrane. The designation Cilium assembly protein DZIP1 Antibody refers to the same gene product.

One major pathway organized through primary cilia is Hedgehog signaling, which contributes to embryonic patterning and tissue maintenance. Experimental work has connected DZIP1 with Hedgehog pathway regulation, although specific interactions can depend on developmental context and organism. The relationship between ciliogenesis and signaling is bidirectional: altered ciliary architecture can change pathway output, while developmental cues influence ciliary dynamics. NSJ Bioreagents includes DZIP1 among proteins involved in centrosome-associated assembly, primary cilium formation, and the spatial organization of developmental signaling. A DZIP1 Antibody can help contextualize the molecular functions and regulation of this protein.

Explore related research targets in our [Developmental Biology Antibodies](#) collection.

Application Notes

The recommended approach is to establish the optimal working dilution of the DZIP1 Antibody / Cilium assembly protein DZIP1 Antibody through preliminary testing.

Immunogen

The DZIP1 Antibody was generated using e. coli-derived human DZIP1 recombinant protein (amino acids R70-D791).

Storage

After reconstitution, keep the DZIP1 Antibody at 4oC for up to one month. Alternatively, store aliquots at -20oC for up to six months. Avoid multiple freezing and thawing cycles.

Alternate Names

DZIP1 antibody, Cilium assembly protein DZIP1 antibody, DAZ interacting zinc finger protein 1 antibody, DZIP antibody, DAZ-interacting zinc finger protein 1 antibody, DAZ-interacting protein 1/2 antibody, DAZ interacting protein 1/2 antibody, DAZ interacting protein 1 antibody