


E2F1 Antibody (R35261)

Catalog No.	Formulation	Size
R35261-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1869
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This E2F1 antibody is available for research use only.



Western blot testing of HeLa nuclear lysate with E2F1 antibody at 1ug/ml. Predicted molecular weight: 48-70 kDa.

Description

Additional name(s) for this target protein: E2F transcription factor 1

Application Notes

Optimal dilution of the E2F1 antibody should be determined by the researcher.

Immunogen

Amino acids QEVTSSEENRATDS were used as the immunogen for this E2F1 antibody.

Storage

Aliquot and store the E2F1 antibody at -20oC.