


E2F6 Antibody (R34878)

Catalog No.	Formulation	Size
R34878-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1876
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This E2F6 antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with E2F6 antibody at 1ug/ml.
Predicted molecular weight: ~32kDa.

Description

Additional name(s) for this target protein: E2F transcription factor 6, E2F-6

Application Notes

Optimal dilution of the E2F6 antibody should be determined by the researcher.

Immunogen

Amino acids NPQQSEELLEVS_N were used as the immunogen for this E2F6 antibody.

Storage

Aliquot and store the E2F6 antibody at -20oC.