


E2F6 Antibody (R34878)

Catalog No.	Formulation	Size
R34878-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1876
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This E2F6 antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with E2F6 antibody at 1ug/ml.
Predicted molecular weight: ~32kDa.

Description

Additional name(s) for this target protein: E2F transcription factor 6, E2F-6

Application Notes

Optimal dilution of the E2F6 antibody should be determined by the researcher.

Immunogen

Amino acids NPQQSEELLEVS_N were used as the immunogen for this E2F6 antibody.

Storage

Aliquot and store the E2F6 antibody at -20oC.