


EDG3 Antibody (R36075)

Catalog No.	Formulation	Size
R36075-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1903
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This EDG3 antibody is available for research use only.



IHC testing of FFPE human kidney with EDG3 antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Endothelial differentiation gene-3; sphingosine-1-phosphate receptor 3 (S1PR3)

Application Notes

Optimal dilution of the EDG3 antibody should be determined by the researcher.

Immunogen

Amino acids SSRKVANHNNSEERS were used as the immunogen for this EDG3 antibody.

Storage

Aliquot and store the EDG3 antibody at -20oC.