


EGR4 Antibody (R34770)

Catalog No.	Formulation	Size
R34770-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1961
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This EGR4 antibody is available for research use only.



Western blot testing of human hippocampus lysate with EGR4 antibody at 0.3ug/ml.
Predicted molecular weight: ~51kDa.

Description

Additional name(s) for this target protein: Early growth response 4

Application Notes

Optimal dilution of the EGR4 antibody should be determined by the researcher.

Immunogen

Amino acids NCGSQGDYQAAPE were used as the immunogen for this EGR4 antibody.

Storage

Aliquot and store the EGR4 antibody at -20oC.