


EIF3E Antibody / eIF-3 p48 (R34522)

Catalog No.	Formulation	Size
R34522-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3646
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This EIF3E antibody is available for research use only.



Western blot testing of Daudi lysate with EIF3E antibody at 0.03ug/ml.
Predicted/observed molecular weight: ~52/48-52kDa.

Description

Additional name(s) for this target protein: Eukaryotic translation initiation factor 3, subunit E, 48kDa

Application Notes

Optimal dilution of the EIF3E antibody should be determined by the researcher.

Immunogen

Amino acids KLNQNSRSEAPN were used as the immunogen for this EIF3E antibody.

Storage

Aliquot and store the EIF3E antibody at -20oC.