


## EIF3M Antibody (R34769)


Catalog No.	Formulation	Size
R34769-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10480
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This EIF3M antibody is available for research use only.



Western blot testing of A549 lysate with EIF3M antibody at 0.1ug/ml. Predicted molecular weight: ~42kDa.



IHC staining of FFPE human testis with EIF3M antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

## Description

Additional name(s) for this target protein: B5 Receptor, PCID1

## Application Notes

Optimal dilution of the EIF3M antibody should be determined by the researcher.

## Immunogen

Amino acids NAWKQNLNKVKN were used as the immunogen for this EIF3M antibody.

## Storage

Aliquot and store the EIF3M antibody at -20oC.