


EIF5A Antibody (R35713)

Catalog No.	Formulation	Size
R35713-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1984
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This eIF5A antibody is available for research use only.



Western blot testing of MOLT4 lysate with eIF5A antibody at 0.1ug/ml. Predicted molecular weight: ~20kDa.

Description

Additional name(s) for this target protein: Eukaryotic translation initiation factor 5A

Application Notes

Optimal dilution of the eIF5A antibody should be determined by the researcher.

Immunogen

Amino acids RPPHRASFLKRLESK were used as the immunogen for this eIF5A antibody.

Storage

Aliquot and store the eIF5A antibody at -20oC.