


## ENDOGL1 Antibody / Endonuclease G-like 1 / EXOG [clone 28E53] (RQ8401)

Catalog No.	Formulation	Size
RQ8401	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Recombinant Rabbit Monoclonal
<b>Isotype</b>	Rabbit IgG
<b>Clone Name</b>	28E53
<b>Purity</b>	Affinity chromatography
<b>UniProt</b>	Q9Y2C4
<b>Applications</b>	Western Blot : 1:500-1:2000
<b>Limitations</b>	This ENDOGL1 antibody is available for research use only.



Western blot testing of human HepG2 cell lysate with ENDOGL1 antibody. Predicted molecular weight ~41 kDa.

### Description

Endo/exonuclease with nicking activity towards supercoiled DNA, a preference for single-stranded DNA and 5'-3' exonuclease activity. [UniProt]

### Application Notes

Optimal dilution of the ENDOGL1 antibody should be determined by the researcher.

## **Immunogen**

A synthetic peptide derived from human ENDOGL1 protein was used as the immunogen for the ENDOGL1 antibody.

## **Storage**

Store the ENDOGL1 antibody at -20°C.