


ENOX2 Antibody (R36443)

Catalog No.	Formulation	Size
R36443-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10495
Localization	Cytoplasmic, membranous
Applications	IHC (FFPE) : 2-4ug/ml
Limitations	This ENOX2 antibody is available for research use only.



IHC testing of FFPE human kidney with ENOX2 antibody at 4ug/ml. HIER: steamed with pH9 Tris/EDTA buffer, HRP-staining.

Description

Additional name(s) for this target protein: Ecto NOX disulfide thiol exchanger 2; tNOX; COVA1; Cytosolic ovarian carcinoma antigen 1; APK1

Application Notes

Optimal dilution of the ENOX2 antibody should be determined by the researcher.

Immunogen

Amino acids EKLKDDKLQVEK were used as the immunogen for this ENOX2 antibody.

Storage

Aliquot and store the ENOX2 antibody at -20oC.