


EPB41L2 Antibody (R33678)

Catalog No.	Formulation	Size
R33678-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2037
Applications	Immunohistochemistry (FFPE) : 2-5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This EPB41L2 antibody is available for research use only.



IHC testing of FFPE human uterus with EPB41L2 at 2ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Erythrocyte membrane protein band 4.1-like 2

Application Notes

Optimal dilution of the EPB41L2 antibody should be determined by the researcher.

Immunogen

Amino acids DYDPEEHGSID were used as the immunogen for this EPB41L2 antibody.

Storage

Aliquot and store the EPB41L2 antibody at -20oC.