

EPB41L3 Antibody (R34214)

Catalog No.	Formulation	Size
R34214-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23136
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This EPB41L3 antibody is available for research use only.



Western blot testing of human brain (cerebellum) lysate with EPB41L3 antibody at 0.3ug/ml. Predicted molecular weight: ~121 kDa, observed here at ~150 kDa.

Description

Additional name(s) for this target protein: Erythrocyte membrane protein band 4.1-like 3, DAL1

Application Notes

Optimal dilution of the EPB41L3 antibody should be determined by the researcher.

Immunogen

Amino acids HKETEITPEDGED were used as the immunogen for this EPB41L3 antibody.

Storage

Aliquot and store the EPB41L3 antibody at -20oC.