


EPB41L3 Antibody (R34214)

Catalog No.	Formulation	Size
R34214-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23136
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This EPB41L3 antibody is available for research use only.



Description

Additional name(s) for this target protein: Erythrocyte membrane protein band 4.1-like 3, DAL1

Application Notes

Optimal dilution of the EPB41L3 antibody should be determined by the researcher.

Immunogen

Amino acids HKETEITPEDGED were used as the immunogen for this EPB41L3 antibody. **Storage** Aliquot and store the EPB41L3 antibody at -20oC.