


## Ephrin A1 Antibody / EFNA1 (R34693)

Catalog No.	Formulation	Size
R34693-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P20827
Gene ID	1942
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Ephrin A1 antibody is available for research use only.



Western blot testing of mouse thymus lysate with Ephrin A1 antibody at 0.01ug/ml.  
Predicted molecular weight: ~24kDa.

## Description

Additional name(s) for this target protein: EPH-related receptor tyrosine kinase ligand 1; LERK-1; Immediate early response protein B61; Tumor necrosis factor alpha-induced protein 4; TNF alpha-induced protein 4

## Application Notes

Optimal dilution of the Ephrin A1 antibody should be determined by the researcher.

1. This Ephrin A1 antibody is specific for isoform a.

## Immunogen

Amino acids KVTVSGKITHSP were used as the immunogen for this Ephrin A1 antibody.

## Storage

Aliquot and store the Ephrin A1 antibody at -20oC.