


ERG Antibody / Transcriptional regulator ERG (R35088)

Catalog No.	Formulation	Size
R35088-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2078
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ERG antibody is available for research use only.



ERG Antibody MOLT4 WB. Western blot testing of human MOLT4 cell lysate with ERG antibody at 0.1ug/ml. Predicted molecular weight ~54 kDa.

Description

Additional name(s) for this target protein: V-ets avian erythroblastosis virus E26 oncogene related; ETS transcription factor ERG

A [microarray-validated ERG antibody](#) for IHC applications is available.

Application Notes

Optimal dilution of the ERG antibody should be determined by the researcher.

Immunogen

Amino acids PNTRLPTSHMP SH were used as the immunogen for this ERG antibody.

Storage

Aliquot and store the ERG antibody at -20oC.