


ERp57 Antibody / PDIA3 / Protein disulfide-isomerase A3 (R33750)

Catalog No.	Formulation	Size
R33750-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2923
Localization	Cytoplasmic, membrane
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This ERp57 antibody is available for research use only.



Western blot testing of HepG2 lysate with ERp57 antibody at 1ug/ml. Predicted molecular weight: ~57-60 kDa.

Description

Additional name(s) for this target protein: GRP58; PDIA3

Application Notes

Optimal dilution of the ERp57 antibody should be determined by the researcher.

Immunogen

Amino acids EEKPKKKKKKAQEDL were used as the immunogen for this ERp57 antibody.

Storage

Aliquot and store the ERp57 antibody at -20oC.