


## ERV3 Antibody (R36494)

Catalog No.	Formulation	Size
R36494-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q14264
Localization	Secreted
Applications	IHC (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This ERV3 antibody is available for research use only.



IHC testing of FFPE human placental tissue with ERV3 antibody at 3.75ug/ml. HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

### Description

Alternate name(s) for this target: Endogenous retroviral sequence 3; Endogenous retrovirus group 3 member 1 Env polypeptide precursor

### Application Notes

Optimal dilution of the ERV3 antibody should be determined by the researcher.

## **Immunogen**

Amino acids PQDNFTLTASS were used as the immunogen for this ERV3 antibody.

## **Storage**

Aliquot and store the ERV3 antibody at -20oC.