

ESR1 Antibody / Estrogen receptor alpha [clone MSVA-564R] (V6071)

Catalog No.	Formulation	Size
V6071-100UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V6071-20UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	MSVA-564R
UniProt	P03372
Localization	Nucleus
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This ESR1/Estrogen receptor alpha antibody is available for research use only.

Application Notes

1. Optimal dilution of the ESR1/Estrogen receptor alpha antibody should be determined by the researcher.
2. This ESR1/Estrogen receptor alpha antibody is recombinantly produced by expression in human HEK293 cells.

Immunogen

A recombinant fragment (around amino acids 495-595) of human ER alpha protein (exact sequence is proprietary) was used as the immunogen for the ESR1/Estrogen receptor alpha antibody.

Storage

ESR1/Estrogen receptor alpha antibody with sodium azide - store at 2 to 8oC; antibody without sodium azide - store at -20 to -80oC.

