

ESR2 Antibody / ER beta [clone ESR2/9710] (V5882)

Catalog No.	Formulation	Size
V5882-100UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V5882-20UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug
V5882SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, kappa
Clone Name	ESR2/9710
UniProt	Q92731
Localization	Nucleus
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This ESR2/ER beta antibody is available for research use only.

Application Notes

Optimal dilution of the ESR2/ER beta antibody should be determined by the researcher.

Immunogen

A recombinant fragment (around amino acids 1-200) of human ESR2 protein (exact sequence is proprietary) was used as the immunogen for the ESR2/ER beta antibody.

Storage

ESR2/ER beta antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

