

ETFB Antibody Biotin Conjugate (R35125BTN)

Catalog No.	Formulation	Size
R35125BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2109
Applications	Western Blot : 0.03-0.1ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:8000
Limitations	This ETFB antibody is available for research use only.



Western blot testing of mouse liver lysate with biotinylated ETFB antibody at 0.03ug/ml.
Predicted molecular weight: ~28kDa.

Description

Additional name(s) for this target protein: Electron-transfer-flavoprotein beta

Application Notes

Optimal dilution of the ETFB antibody should be determined by the researcher.

Immunogen

Amino acids PQRTAGVKVETTE were used as the immunogen for this ETFB antibody.

Storage

Aliquot and store the ETFB antibody at -20oC.