


## Eukaryotic Translation Initiation Factor 6 / EIF6 (R32150)


| Catalog No. | Formulation                                           | Size   |
|-------------|-------------------------------------------------------|--------|
| R32150      | 0.5mg/ml if reconstituted with 0.2ml sterile DI water | 100 ug |

**Bulk quote request**


|                           |                                                                                              |
|---------------------------|----------------------------------------------------------------------------------------------|
| <b>Availability</b>       | 1-3 business days                                                                            |
| <b>Species Reactivity</b> | Human, Mouse, Rat                                                                            |
| <b>Format</b>             | Antigen affinity purified                                                                    |
| <b>Host</b>               | Rabbit                                                                                       |
| <b>Clonality</b>          | Polyclonal (rabbit origin)                                                                   |
| <b>Isotype</b>            | Rabbit IgG                                                                                   |
| <b>Purity</b>             | Antigen affinity                                                                             |
| <b>Buffer</b>             | Lyophilized from 1X PBS with 2.5% BSA, 0.025% sodium azide                                   |
| <b>UniProt</b>            | P56537                                                                                       |
| <b>Localization</b>       | Nuclear and cytoplasmic                                                                      |
| <b>Applications</b>       | Western Blot : 0.1-0.5ug/ml<br>IHC (FFPE) : 0.5-1ug/ml<br>ICC (FFPE) : 0.5-1ug/ml            |
| <b>Limitations</b>        | This Eukaryotic Translation Initiation Factor 6 antibody is available for research use only. |




Western blot testing of 1) rat heart, 2) rat liver, 3) mouse liver, 4) human placenta, 5) COLO320 and 6) HeLa lysate with Eukaryotic Translation Initiation Factor 6 antibody. Expected molecular weight ~27 kDa.




IHC testing of FFPE human intestinal cancer tissue with Eukaryotic Translation Initiation Factor 6 antibody. HIER: Boil the paraffin sections in pH 6, 10mM citrate buffer for 20 minutes and allow to cool prior to staining.




IHC testing of FFPE mouse intestine with Eukaryotic Translation Initiation Factor 6 antibody. HIER: Boil the paraffin sections in pH 6, 10mM citrate buffer for 20 minutes and allow to cool prior to staining.



IHC testing of FFPE rat testis with Eukaryotic Translation Initiation Factor 6 antibody. HIER: Boil the paraffin sections in pH 6, 10mM citrate buffer for 20 minutes and allow to cool prior to staining.



IHC testing of FFPE human SMMC-7221 cells with Eukaryotic Translation Initiation Factor 6 antibody. HIER: Boil the paraffin sections in pH 6, 10mM citrate buffer for 20 minutes and allow to cool prior to staining.



IHC testing of FFPE human SW480 cells with Eukaryotic Translation Initiation Factor 6 antibody. HIER: Boil the paraffin sections in pH 6, 10mM citrate buffer for 20 minutes and allow to cool prior to staining.

## Description

EIF6 (Eukaryotic Translation Initiation Factor 6), also called EIF3A or ITGB4BP, is a human gene. By fluorescence in situ hybridization, Sanvito et al. (1998) mapped the ITGB4BP gene to 20q11.2. Ceci et al. (2003) demonstrated that the ribosomal 60S subunit is activated by release of EIF6. In the cytoplasm, EIF6 is bound to free 60S but not to 80S subunits. Furthermore, EIF6 interacts in the cytoplasm with RACK1, a receptor for activated protein kinase C. Gandin et al. (2008) demonstrated that mammalian eIF6 is required for efficient initiation of translation in vivo. Eif6-null mouse embryos were lethal at preimplantation. Heterozygous mice had 50% reduction of eIF6 levels in all tissues, and showed reduced mass of hepatic and adipose tissues due to a lower number of cells and to impaired G1/S cell cycle progression.

## Application Notes

Optimal dilution of the Eukaryotic Translation Initiation Factor 6 antibody should be determined by the researcher.

## Immunogen

Amino acids 66-210 of the human protein were used as the immunogen for the Eukaryotic Translation Initiation Factor 6 antibody.

## Storage

After reconstitution, the Eukaryotic Translation Initiation Factor 6 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.