


EVC2 Antibody (R34729)

Catalog No.	Formulation	Size
R34729-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	132884
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This EVC2 antibody is available for research use only.



IHC testing of FFPE human placenta with EVC2 antibody at 2 ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Ellis van Creveld syndrome 2, LBN, LIMBIN

Application Notes

Optimal dilution of the EVC2 antibody should be determined by the researcher.

1. This EVC2 antibody will detect isoforms 1 and 2.

Immunogen

Amino acids LNAKKAMRALGMD were used as the immunogen for this EVC2 antibody.

Storage

Aliquot and store the EVC2 antibody at -20oC.