


EXO1 Antibody (R34271)

Catalog No.	Formulation	Size
R34271-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9156
Applications	Western Blot : 0.2-0.5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This EXO1 antibody is available for research use only.



Western blot testing of A549 cell lysate with EXO1 antibody at 0.2ug/ml. The expected 94-115 kDa band and the additional ~38 kDa band are both blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: Exonuclease 1

Application Notes

Optimal dilution of the EXO1 antibody should be determined by the researcher.

Immunogen

Amino acids HRNYSPRPESGT were used as the immunogen for this EXO1 antibody.

Storage

Aliquot and store the EXO1 antibody at -20oC.