


EYA1 Antibody (R36145)

Catalog No.	Formulation	Size
R36145-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2138
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This EYA1 antibody is available for research use only.



Western blot testing of HEK293 lysate with EYA1 antibody at 2ug/ml. Predicted molecular weight: ~61kDa.

Description

Additional name(s) for this target protein: Eyes absent homolog 1

Application Notes

Optimal dilution of the EYA1 antibody should be determined by the researcher.

Immunogen

Amino acids TDPTAEYSTIHSP were used as the immunogen for this EYA1 antibody.

Storage

Aliquot and store the EYA1 antibody at -20oC.