


## FABP3 Antibody (cardiac) (R35172)

| Catalog No.  | Formulation   | Size   |
|--------------|---|--------|
| R35172-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |   |
|----------------------|---|
| Availability         | 1-3 business days   |
| Species Reactivity   | Human, Mouse, Rat   |
| Predicted Reactivity | Cow   |
| Format               | Antigen affinity purified                                       |
| Clonality            | Polyclonal (goat origin)  |
| Isotype              | Goat Ig   |
| Purity               | Antigen affinity  |
| Gene ID              | 2170  |
| Localization         | Cytoplasmic, membrane   |
| Applications         | Western Blot : 0.01-0.03ug/ml<br>ELISA (peptide) LOD : 1:128000 |
| Limitations          | This FABP3 antibody is available for research use only.         |



Western blot testing of human [A], mouse [B] and rat [C] heart lysates with FABP3 antibody at 0.01ug/ml. Predicted molecular weight: ~15 kDa.

### Description

Additional name(s) for this target protein: Fatty acid binding protein 3, heart-type FABP

### Application Notes

Optimal dilution of the FABP3 antibody should be determined by the researcher.

## Immunogen

Amino acids PTIIIEKNGDILTLK were used as the immunogen for this FABP3 antibody.

## Storage

Aliquot and store the FABP3 antibody at -20oC.