


FABP3 Antibody (cardiac) (R35172)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35172-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|-----------------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human, Mouse, Rat |
| Predicted Reactivity | Cow |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 2170 |
| Localization | Cytoplasmic, membrane |
| Applications | Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:128000 |
| Limitations | This FABP3 antibody is available for research use only. |



Western blot testing of human [A], mouse [B] and rat [C] heart lysates with FABP3 antibody at 0.01ug/ml. Predicted molecular weight: ~15 kDa.

Description

Additional name(s) for this target protein: Fatty acid binding protein 3, heart-type FABP

Application Notes

Optimal dilution of the FABP3 antibody should be determined by the researcher.

Immunogen

Amino acids PTTIIEKNGDILTLK were used as the immunogen for this FABP3 antibody.

Storage

Aliquot and store the FABP3 antibody at -20°C.