


FABP4 Antibody (R36099)

Catalog No.	Formulation	Size
R36099-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2167
Localization	Cytoplasmic
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This FABP4 antibody is available for research use only.



FABP4 Antibody Adipose Tissue WB. Western blot testing of human adipose lysate with FABP4 antibody at 0.03ug/ml. Predicted molecular weight ~15 kDa.

Description

Additional name(s) for this target protein: Fatty acid binding protein 4

Additional studies involving adipocyte biology, lipid transport signaling, and metabolic-inflammatory pathway regulation

may benefit from our [FABP4 antibody](#) page featuring clone FABP4/4424 with HuProt protein microarray specificity validation.

Application Notes

Optimal dilution of the FABP4 antibody should be determined by the researcher.

Immunogen

Amino acids SVNGDVITIKSEST were used as the immunogen for this FABP4 antibody.

Storage

Aliquot and store the FABP4 antibody at -20oC.