


## FABP5 Antibody (epidermal) (RQ4416)


Catalog No.	Formulation	Size
RQ4416	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)


<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit IgG
<b>Purity</b>	Antigen affinity purified
<b>Buffer</b>	Lyophilized from 1X PBS with 2% Trehalose
<b>UniProt</b>	Q01469
<b>Localization</b>	Cytoplasmic, nuclear
<b>Applications</b>	Western Blot : 0.5-1ug/ml Immunofluorescence : 5ug/ml
<b>Limitations</b>	This FABP5 antibody is available for research use only.




FABP5 Antibody U-2 OS Cell IF. Immunofluorescent staining of FFPE human U-2 OS cells with FABP5 antibody (green) and DAPI nuclear stain (blue). HIER: steam section in pH6 citrate buffer for 20 min.



FABP5 Antibody A549 Cell IF. Immunofluorescent staining of FFPE human A549 cells with FABP5 antibody (green) and DAPI nuclear stain (blue). HIER: steam section in pH6 citrate buffer for 20 min.



Western blot testing of 1) human HeLa, 2) human A549, 3) rat thymus and 4) mouse thymus lysate with FABP5 antibody at 0.5ug/ml. Predicted molecular weight ~15 kDa.



Western blot testing of 1) human RT4, 2) human PC-3, 3) rat lung and 4) mouse lung tissue lysate with FABP5 antibody at 0.5ug/ml. Predicted molecular weight ~15 kDa.

## Description

FABP5 Antibody detects Fatty acid-binding protein, epidermal, a protein that in humans is encoded by the FABP5 gene. This gene encodes the fatty acid binding protein found in epidermal cells, and was first identified as being upregulated in psoriasis tissue. It is mapped to 8q21.13. Fatty acid binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. It is thought that FABPs roles include fatty acid uptake, transport, and metabolism.

Explore our central [FABP5 antibody](#) page for additional western blot, immunohistochemistry, and microarray specificity validation data supporting studies of lipid transport, epidermal biology, and metabolic signaling pathways.

## Application Notes

Optimal dilution of the FABP5 antibody should be determined by the researcher.

## Immunogen

Amino acids KWRLMESHGFEEYMKELGVGLALRKMAAMAKPD from the human protein were used as the immunogen for the FABP5 antibody.

## Storage

After reconstitution, the FABP5 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.